

---

## ***Biomphalaria glabrata* embryonic cell (Bge) culture protocol**

*Authors: Wannaporn Ittiprasert, PhD  
André Miller*

### **Materials**

Bge Medium 1 liter:

- Schneiders's Drosophila Medium 220 mL
- lactalbumin hydrolysate 4.5 g
- Galactose 1.3g
- MilliQ H2O up to 900mL
- Gentamycin (10 mg/mL stock) 2.0 mL
- Phenol Red (0.5% solution) 1.62mL

Adjust pH to 7.0

Sterile filter at least > 0.22 micron

### **Grow cells at 26°C in complete Bge medium (7-10% heat-inactivated FCS in Bge medium)**

Keep the cells growing at a high density. This cell line likes to form cell balls on top of a monolayer. Split the cells when the cell balls get very dense or about once a week. Split the cells only 1:2, keeping the cell culture dense.

**To split/passage Bge cells:** remove the attached cells by using a cell scraper (not trypsin) and transfer the suspension into a 15 ml tube. Resuspend cells in 1 mL complete Bge medium and use 0.1-0.2mL of cell suspension to re-seed in to new 50 mL tissue culture flask with 5 mL complete Bge medium.

**To freeze Bge cells:** remove the attached cells by using a cell scraper and transfer the cell suspension into a 15 mL tube. Centrifuge at 700 rpm for 5 min. Re-suspend cells in 0.5-1 mL of complete Bge medium, then count cells by cell counter. Dilute the cells to a final concentration at  $10^9$ - $10^{10}$  cells/mL in room temperature-freezing medium [9 FBS: 1 DMSO (v/v)]. Aliquot 0.8-1.5 mL into freezing vials then transfer into freezing container (pre-cool down in refrigerator; 0-4°C) as soon as possible. Keep freezing container in -70°C for overnight, then move the vial to liquid nitrogen storage.

**To revive frozen cells from liquid nitrogen:** thaw cells in 35-37°C water bath until ~80% is thawed, then add 1-2 ml of Bge medium at room temperature immediately. Transfer the cell suspension into a 15 mL tube and centrifuge at 700 rpm for 5 min, room temperature (RT). Resuspend the cell pellet with at least 7 mL of Bge medium and centrifuge as above. Resuspend the cells with 5 mL complete Bge medium, and then transfer the cells to 50 mL-tissue culture flask. Let the cells rest in culture at 26°C for a week (DO NOT CHANGE MEDIUM) to let cells adapt and start dividing before changing medium.

### **Comments**

Bge cells are not easily retrieved from frozen vials, therefore it is best to keep a culture growing at all times. For more information, contact André Miller at [amiller@afbr-bri.com](mailto:amiller@afbr-bri.com)